



# UNLOCK THE ANSWER

The Reveal LINQ™ System  
lets your physician learn about  
your heart while you live your life.



Medtronic

# WHY LONG-TERM MONITORING MAY BE RIGHT FOR YOU

**Accurate information is the key to better healthcare.**

But irregular heartbeats may occur infrequently, making the underlying condition hard to find and treat.

**Now, there's a way.** With continuous, long-term monitoring, your doctor can get information about your heart to help find the right answer for you.



## **UNEXPLAINED FAINTING—**

a sudden loss of consciousness—also called “passing out” or syncope—that occurs when blood pressure drops and oxygen can’t reach the brain.

- **Why monitor?** Fainting can be serious when caused by an abnormal heart rhythm.

## **HEART PALPITATIONS—**

a usually harmless condition in which the heart feels as if it’s pounding, racing or fluttering.

- **Why monitor?** Palpitations may indicate a serious underlying heart condition.

## **ATRIAL FIBRILLATION—**

a common condition in which the upper chambers of the heart beat very fast and irregularly. Also called AF or AFib.

- **Why monitor?** AFib requires treatment because it’s strongly linked to stroke.

## **UNEXPLAINED STROKE—**

a stroke of unknown cause that may be triggered by AFib, causing a blood clot to form and block oxygen to the brain.

- **Why monitor?** Detecting AFib after stroke may help physicians prevent a second one.

# THE REVEAL LINQ™ SYSTEM

The Reveal LINQ System includes a small Insectable Cardiac Monitor (or ICM) placed just under the skin during a minimally-invasive outpatient procedure; and the MyCareLink™ Patient Monitor—a bedside unit that collects heart rhythm data from the ICM and makes it available to your doctor.

- **Discreet**—The Reveal LINQ ICM is not visible in most patients.
- **Easy**—The MyCareLink Monitor uses a cellular signal and a power outlet.
- **Portable**—Go wherever you want. You can travel and take the MyCareLink Monitor with you.
- **MRI Compatible**—The Reveal LINQ ICM is safe for use in an MRI setting.
- **3 Years**—The Reveal LINQ ICM detects abnormal heart rhythms for up to 3 years.

# MARK SYMPTOMS AS THEY HAPPEN

Your doctor may choose to add a Patient Assistant to your Reveal LINQ System if he or she wants to compare specific heart rhythm data to symptoms you experience.

The Patient Assistant is a small device used to mark symptoms in the Reveal LINQ ICM. The MyCareLink Monitor will then send this information during your next automatic transmission.

Medtronic is committed to partnering with your doctor and unlocking the answers to your health.



Reveal LINQ ICM



MyCareLink  
Patient Monitor



Patient Assistant

# HOW THE REVEAL LINQ SYSTEM WORKS

Once your Reveal LINQ ICM is inserted, all you need to do is leave the MyCareLink Monitor plugged in, powered on and in a location that receives an adequate cellular signal.

**THE SYSTEM** will take care of the rest!



The process for sending information is easy and worry-free.



Heart data is recorded with the Reveal LINQ ICM.



You go to bed near your MyCareLink Monitor.



Your data is shared with the MyCareLink Monitor.



Your doctor can access your data.



Your doctor will contact you if necessary.


# THINGS TO CONSIDER

All surgical procedures carry risk, but those associated with the Reveal LINQ System are rare. Since the ICM is inserted just under the skin, there is a small risk of infection and/or sensitivity. **Talk to your doctor about any risks particular to your health and situation.**

## WHY REVEAL LINQ?

There are many benefits to the Reveal LINQ System:

- Insertion procedure is minimally-invasive, easy and brief
- Long-term monitoring works continuously
- Ability to capture irregular heartbeats that are infrequent or unpredictable helps with diagnosis
- Use of Reveal LINQ System doesn't require a change in daily activities
- Reveal LINQ ICM is safe for use in an MRI setting



# GET MORE INFORMATION

Medtronic is dedicated to providing you with the most up-to-date information about the Reveal LINQ System.



Learn more about long-term heart monitoring, read patient stories and find answers about the Reveal LINQ System at:

[MonitorYourHeart.com/LINQ](http://MonitorYourHeart.com/LINQ)

Patient service specialists are available by phone from 7 a.m. to 7 p.m. (CST), Monday through Friday.

**1 (800) 551-5544**



Please talk to your doctor if you are interested in finding out more about the Reveal LINQ System or if you are ready to schedule your insertion.





# INFORMATION SENT AUTOMATICALLY TO YOU

The Reveal LINQ System educational emails are an easy way for you to stay informed about long-term heart monitoring with the Reveal LINQ System.

## Sign up to receive emails that will help you:

- Learn more about the purpose of long-term heart monitoring
- Set up and use your Reveal LINQ System
- Understand the system components and what they do
- Receive relevant FAQs, answering questions as they come up

**SIGN UP** by returning this mail-in form  
**OR** fill out our online form at  
[MonitorYourHeart.com/LINQsignup](https://www.monitoryourheart.com/LINQsignup)

---

Name

---

Email

---

Signature

---

Scheduled date for Reveal LINQ ICM insertion


Please select the reason you are considering the Reveal LINQ System

- |   |  |
|---|--|
| <input type="checkbox"/> Unexplained fainting | <input type="checkbox"/> Heart palpitations            |
| <input type="checkbox"/> Atrial fibrillation  | <input type="checkbox"/> Suspected atrial fibrillation |
| <input type="checkbox"/> Unexplained stroke   | <input type="checkbox"/> Not sure                      |


By completing and submitting this form, you are granting Medtronic permission to add the personal information you provide on this form to Medtronic's patient database, which Medtronic will use to communicate directly with you. You also agree to receive non-password protected emails from Medtronic for the purpose of providing educational information in support of product use or to inquire about your experience regarding our products or services.

Medtronic respects the confidentiality of your personal information. We will not share your personal information except as described above. If, at any time, you wish to revoke all or part of this permission, you can follow the "Unsubscribe" link in any email you receive from us or send a request in writing to:

Reveal LINQ Educational Emails  
912 Industrial Drive  
Milbank, SD 57252




NO POSTAGE  
NECESSARY  
IF MAILED  
IN THE  
UNITED STATES



**BUSINESS REPLY MAIL**  
FIRST-CLASS MAIL PERMIT NO. 7611 MINNEAPOLIS, MN

POSTAGE WILL BE PAID BY ADDRESSEE

REVEAL LINQ EDUCATIONAL EMAILS  
912 INDUSTRIAL DR  
MILBANK SD 57252-9920





## Medtronic CareLink™ Monitor and Medtronic CareLink Network

**Intended Use:** The CareLink™ Monitor and the CareLink Network are indicated for use in the transfer of patient data from some Medtronic implantable cardiac devices based on physician instructions and as described in the product manual. These products are not a substitute for appropriate medical attention in the event of an emergency and should only be used as directed by a physician.

**Contraindications:** There are no contraindications for the CareLink Monitor.

**Warnings and Precautions:** The CareLink Monitor must only be used for interrogating compatible Medtronic implantable devices. The CareLink Monitor is intended for use within the prescribing country.

## Reveal LINQ™ LNQ11 Insertable Cardiac Monitor and Patient Assistant

The Reveal LINQ Insertable Cardiac Monitor is an implantable patient-activated and automatically-activated monitoring system that records subcutaneous ECG and is indicated in the following cases:

- patients with clinical syndromes or situations at increased risk of cardiac arrhythmias
- patients who experience transient symptoms such as dizziness, palpitation, syncope, and chest pain, that may suggest a cardiac arrhythmia

The device has not been tested specifically for pediatric use.

Possible risks associated with the implant of the Reveal LINQ Insertable Cardiac Monitor include, but are not limited to, infection at the surgical site, device migration, erosion of the device through the skin and/or sensitivity to the device material.

Treatment with a Reveal LINQ Insertable Cardiac Monitor is prescribed by your physician. This treatment is not for everyone. Please talk to your doctor to see if it is right for you. Your physician should discuss all potential benefits and risks with you. Although many patients benefit from the use of this treatment, results may vary.

*See the device manual for detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential complications/adverse events. For further information, please call Medtronic at 1 (800) 328-2518 and/or consult Medtronic's website at [www.medtronic.com](http://www.medtronic.com).*

**Caution:** Federal law (USA) restricts these devices to sale by or on the order of a physician.

# Medtronic

### Patient Services

Medtronic  
8200 Coral Sea St. NE MVS14  
Mounds View, MN 55112  
Patient toll-free line:  
1 (800) 551-5544  
Fax: (763) 367-5809  
7:00 am to 7:00 pm CT Mon.-Fri.

UC201602395 EN © Medtronic 2015.  
Minneapolis, MN. All Rights Reserved.  
Printed in USA. 10/2015

[medtronic.com](http://medtronic.com)